

Stance - Control Knee™

Designed for people who have knee instability due to quadriceps weakness or absence. Results in a more normal and energy-efficient gait than traditional locked knee joints

SCOKJ® Facts



Prevents knee collapse regardless of flexion


Stance control can be triggered by weightbearing or joint motion, according to patient needs

Provides three operating modes: Automatic Stance Control, Free Motion, and Manual Lock

Suitable for metal, thermoplastic, laminated, and hybrid orthoses

Fabrication and adjustments are straightforward and require no expensive special tooling

Can be used with bilateral knee paralysis or paresis



Polio patient wearing bi-lateral SCOKJs®



Spinal Cord Injury

Horton Technology Inc.

Innovations In Orthotics and Prosthetics

5220 West 12th * Little Rock, AR * 72204

1-866-663-3970

www.stancecontrol.com

Visit us in New Orleans booth # 32